

Kompatibilitetsoversigt for Oticon opladere



Oticon Charger
miniRITE



Oticon SmartCharger
miniRITE R



Oticon SmartCharger
miniBTE R



Charger 1.0
miniRITE R



Charger 1.0
miniBTE R

Oticon Intent™	•	-	-	-	-
Oticon Real™	-	•	•	•	•
Oticon More™	-	•	•	•	•
Oticon Zircon	-	•	•	•	•
Oticon Opn S™	-	-	-	•	-
Oticon Ruby	-	-	-	•	-
Oticon Play PX	-	•	•	•	•
Oticon Opn Play™	-	-	-	•	-
Oticon CROS PX	-	•	-	•	-