

Bernafon Alpha

Better hearing starts now



Bernafon

THE SOUND EXPERIENCE

Meet Bernafon Alpha, your new partner in better hearing



Versatile and practical

Bernafon Alpha comes with many functionalities for easy handling.

2.4 GHz Bluetooth® Low-Energy wireless technology streams sound from external audio devices

LED light indicates power status

Double push buttons change program & volume independently

Wide range of domes and custom molds ensure a secure and comfortable fit

Powered by zinc-air batteries

Dust and water-resistant coating extends durability (IP68-rated)

Telecoil connects to available induction loop systems in public places



Elegant in design

Choose between **5 attractive colors** to match your personal style.



Antique Bronze

Metallic Anthracite

Base Shell: Metallic Anthracite



Metallic Silver



Metallic Anthracite

Sand Beige

Base Shell: Metallic Silver

The Bernafon App

Monitor your hearing aids with the free **Bernafon App**. Use it to change your hearing aids volume or program and check their battery levels. You can also use the "Find my hearing aid" function, should you need it.



Scan me!



Scan me!



App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Accessories for additional support




Clip **SoundClip-A** to your jacket and connect it to your smartphone for hands-free phone calls. Or use it as a remote microphone to improve communication over distances of up to approximately 66 feet.



The **TV-A** lets you stream TV entertainment into your hearing aids. Enjoy the latest action movie in stereo sound quality. Meanwhile, your friends or family can watch TV with you at the volume best suited for them.

Stream sound directly to your ears


Not only a hearing aid, but a technological gadget too. Bernafon Alpha connects to many external audio devices via 2.4 GHz Bluetooth® Low-Energy wireless technology. Answer phone calls, enjoy music, or listen to audio books directly streamed into your hearing aids. They function as high-quality stereo headphones and stream sound for many hours.



For information on compatibility, please visit [bernafon.us/support/compatibility-guide](https://www.bernafon.us/support/compatibility-guide)

Apple, the Apple logo, iPhone, iPad, Mac and the Mac logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Demant A/S is under license. Other trademarks and trade names are those of their respective owners.





The Bernafon Sound Experience

At Bernafon, we make it our job to empower the work of hearing care professionals. That's because professional care is the key that turns extraordinary technology into personal solutions. Solutions that help people with hearing loss live life the way they want to.

As part of the Demant Group, with over a century of experience and know-how, Bernafon benefits from one of the world's strongest research and development teams within hearing healthcare. We bring this expertise and dedication to hearing care professionals – so they can provide an exceptional sound experience.

For more information on our hearing aids, visit our Bernafon website at: [bernafon.us](https://www.bernafon.us).