

# Översikt kompatibilitet laddare Bernafon hörapparater



Bernafon Charger Plus miniRITE      Bernafon Charger miniRITE      Charger Plus miniRITE T R      Charger Plus miniBTE T R      Charger miniRITE T R      Charger miniBTE T R

Bernafon Encanta	•	•	–	–	–	–
Bernafon Alpha XT	–	–	•	•	•	•
Bernafon Alpha	–	–	•	•	•	•
Bernafon Viron	–	–	–	–	•	–
Bernafon CROS		–	•	–	•	–