

Bernafon Alpha XT

A natural sound experience



Bernafon
THE SOUND EXPERIENCE



Embrace the experience of Bernafon Alpha XT.

Available in both rechargeable and non-rechargeable styles, Alpha XT combines exceptional sound quality with advanced connectivity, offering versatility and elegance.

Bernafon Alpha XT hearing aids:

- Simultaneously augments speech understanding and listening comfort
- Lasts all day long*. Rechargeable model charges in under 3 hours.
- Sits comfortably behind the ear for a discreet and secure fit all day long
- Streams sound directly from your smartphone or other audio devices
- Connects to the Bernafon App so you can easily monitor and control your hearing aids and benefit from remote hearing aid fittings


Choose between 10 attractive color combinations to match your personal style



*Lithium-ion battery performance varies depending upon hearing loss, lifestyle and streaming behavior.

Rechargeable for a full day of use*, including sound streaming

With the rechargeable Alpha XT hearing aids, simply put your hearing aids in the charging station overnight and comfortably start your day.



Or make use of the **portable Charger Plus**. Equipped with a power bank, it offers 3 full charging cycles of two hearing aids without the need for a power supply. Recharge the power bank of the Charger Plus in only 4 hours and travel comfortably and hassle-free.



Accessories for additional support



Clip **SoundClip-A** to your jacket and connect it to your smartphone for hands-free phone calls. Or use it as a remote microphone to improve communication over distances of up to approximately 66 feet.



The **TV-A** lets you stream TV entertainment into your hearing aids. Enjoy the latest action movie in stereo sound quality. Meanwhile, your friends or family can watch TV with you at the volume best suited for them.

*Lithium-ion battery performance varies depending upon hearing loss, lifestyle and streaming behavior.

Direct sound streaming

Bernafon Alpha XT hearing aids connect to external audio devices* via 2.4 GHz Bluetooth® Low Energy wireless technology. Answer phone calls, enjoy music, or listen to audio books – your hearing aids function as high-quality stereo headphones and stream sound for many hours.




Benefit from **hands-free communication** when using an iPhone®* or iPad®*.



*For information on compatibility, please visit bernafon.us/support/compatibility-guide

Download the Bernafon App



The Bernafon App enables you to:

- Easily control hearing aid volume and program
- Check on battery level
- Find misplaced hearing aids
- Use Remote Fitting to connect to your hearing care professional for remote appointments

Apple, the Apple logo, iPhone, iPad, Mac and the Mac logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Demant A/S is under license. Other trademarks and trade names are those of their respective owners.



The Bernafon Sound Experience

At Bernafon, we make it our job to empower the work of hearing care professionals. That's because professional care is the key that turns extraordinary technology into personal solutions. Solutions that help people with hearing loss live life the way they want to.

As part of the Demant Group, with over a century of experience and know-how, Bernafon benefits from one of the world's strongest research and development teams within hearing healthcare. We bring this expertise and dedication to hearing care professionals – so they can provide an exceptional sound experience.

For more information on our hearing aids, visit our Bernafon website at: bernafon.us.